

GOAT ANTI-PGLYRP ANTIBODY

SKU: EB05467



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	TNF superfamily, member 3 (LTB)-like (peptidoglycan recognition protein) TNFSF3L PGRP-S TAG7 PGRP peptidoglycan recognition protein PGLYRP
Accession ID	NP_005082
Immunogen	Peptide with sequence C-YHLIQNWPHYRSP, from the C Terminus of the protein sequence according to NP_005082.
Peptide Sequence	C-YHLIQNWPHYRSP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	8993
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 18-22kDa band seen in MOLT-4 lysate. Recommended for use at 2-4µg/ml.