

GOAT ANTI-PHEMX / TSPAN32 ANTIBODY

SKU: EB05430



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	TSPAN32 PHEMX tetraspanin 32 MGC22455 PHMX TSSC6 pan-hematopoietic expression pan-hematopoietic expression protein tetraspanin tumor-suppressing STF cDNA 6 tumor-suppressing subchromosomal transferable fragment cDNA 6 tumor-suppressing subtransferable candidate 6 ART1 FLJ17158 FLJ97586 MGC22455
Accession ID	NP_620591.3
Blocking Peptide	EBP05430
Immunogen	Peptide with sequence GGLSGCPERGLSD, from the C Terminus of the protein sequence according to NP_620591.3.
Peptide Sequence	GGLSGCPERGLSD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	10077
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx 33-35kDa band seen in HeLa lysate. Recommended for use at 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

