

GOAT ANTI-PHF19 (ISOFORM B) ANTIBODY

SKU: EB12288



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	tudor domain containing 19B polycomb-like protein 3 polycomb-like 3 polycomb like 3 hPCL3 TDRD19B PCL3 MTF2L1 PHD finger protein 19 PHF19
Accession ID	NP_001009936.1
Blocking Peptide	EBP12288
Immunogen	Peptide with sequence C-DQRLPSQSLSSKQKGH, from the internal region (near C Terminus) of the protein sequence according to NP_001009936.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_056466.1; NP_001009936.1).
Peptide Sequence	C-DQRLPSQSLSSKQKGH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	26147
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)