Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PHF19 (ISOFORM B) ANTIBODY

SKU: EB12288

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

tudor domain containing 19B|polycomb-like protein 3|polycomb-like 3|polycomb like **Alias**

3|hPCL3|TDRD19B|PCL3|MTF2L1|PHD finger protein 19|PHF19 **Names**

Accession

NP 001009936.1 ID

Blocking

EBP12288 **Peptide**

Peptide with sequence C-DQRLPSQSLSSKQKGH, from the internal region (near C Terminus) of the protein Immunogen

sequence according to NP 001009936.1.

Product Comments

This antibody is expected to recognize both reported isoforms (NP_056466.1; NP_001009936.1).

Peptide

C-DQRLPSQSLSSKQKGH Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human Species

Human

26147

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:128000.

Limit

Application Pep-ELISA **Type**

> Goat Anti-PHF19 (isoform b) Antibody https://everestbiotech.com/products/goat-anti-phf19-isoform-b-antibody/