

GOAT ANTI-PIK3CG ANTIBODY

SKU: EB05591



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	PI3Kgamma phosphatidylinositol 3-kinase, catalytic, gamma polypeptide phosphatidylinositol 3-kinase catalytic 110-kD gamma phosphoinositide-3-kinase gamma catalytic subunit phosphatidylinositol 3 kinase gamma, p110 gamma 1-phosphatidylinositol 3-kinase PTDINS-3-kinase p110-gamma PI3KgammaPI3-kinase PI3CG PIK3 PI3K phosphoinositide-3-kinase, catalytic, gamma polypeptide PIK3CG
Accession ID	NP_002640.2
Blocking Peptide	EBP05591
Immunogen	Peptide with sequence C-VLGIKQGEKHS, from the C Terminus of the protein sequence according to NP_002640.2.
Peptide Sequence	C-VLGIKQGEKHS
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	5294
Mouse Gene ID	30955
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx 110-120kDa band seen in Human Heart lysate (Predicted MW of 121kDa according to NP_002640). Recommended for use at 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa