



Telephone: (650) 697-3600



GOAT ANTI-PIN1 ANTIBODY

SKU: EB05326

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

peptidyl-prolyl cis/trans isomerase, NIMA-interacting|dod|DOD|protein (peptidyl-prolyl cis/trans isomerase) **Alias**

NIMA-interacting 1|PIN1 **Names**

Accession

NM 006221

Immunogen Peptide with sequence C-TDSGIHIILRTE, from the C Terminus of the protein sequence according to NM_006221.

Peptide

C-TDSGIHIILRTE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse **Species**

Reactive Human

Species

Human

Gene ID

5300

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

Western

Western blot: Approx 20-22kDa bands seen in Human Kidney lysate. Recommended for use at 1-3µg/ml. Blot