

GOAT ANTI-PIN1 ANTIBODY

SKU: EB05326



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	peptidyl-prolyl cis/trans isomerase, NIMA-interacting 1 DOD protein (peptidyl-prolyl cis/trans isomerase) NIMA-interacting 1 PIN1
Accession ID	NM_006221
Immunogen	Peptide with sequence C-TDSGIHILRTE, from the C Terminus of the protein sequence according to NM_006221.
Peptide Sequence	C-TDSGIHILRTE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	5300
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 20-22kDa bands seen in Human Kidney lysate. Recommended for use at 1-3µg/ml.