

GOAT ANTI-PINX1 (AA225-237) ANTIBODY

SKU: EB09726



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PINX1 FLJ20565 LPTL LPTS MGC8850 67-11-3 protein hepatocellular carcinoma-related putative tumor suppressor PIN2-interacting protein 1 Pin2-interacting protein X1
Accession ID	NP_060354.4
Blocking Peptide	EBP09726
Immunogen	Peptide with sequence C-ESYLQPKAKRHTE, from the internal region of the protein sequence according to NP_060354.4.
Peptide Sequence	C-ESYLQPKAKRHTE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	54984
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 40kDa band observed in lysates of cell lines A549, HepG2 and K562 (calculated MW of 37.0kDa according to NP_060354.4). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

