Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PIRIN ANTIBODY

SKU: EB05026

×

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

 $\textbf{Unit Size} \qquad 100 \; \mu \text{g}$

Storage

Instructions

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Pirin|PIR

Names

Accession NP_003653

Peptide with sequence C-FERAKTWKSKIGN, from the C Terminus of the protein sequence according to Immunogen

NP 003653.

Peptide

Sequence C-FERAKTWKSKIGN

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions Refrigerated

Predicted

Species

Human, Mouse

Reactive Species

Human

Human

8544

Gene ID Product

Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Western

Blot Western blot: Approx 35kDa band seen in HeLa cell lysate. Recommended for use at 4-6 μg/ml.