

GOAT ANTI-PIRIN ANTIBODY

SKU: EB05026

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	Pirin PIR
Names	
Accession ID	NP_003653
Immunogen	Peptide with sequence C-FERAKTWKSKIGN, from the C Terminus of the protein sequence according to NP_003653.
Peptide Sequence	C-FERAKTWKSKIGN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	8544
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 35kDa band seen in HeLa cell lysate. Recommended for use at 4-6 µg/ml.