

GOAT ANTI-PITPN / PITP ALPHA ANTIBODY

SKU: EB06319



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names phosphatidylinositol transfer protein|phosphotidylinositol transfer protein|PITPNB|PITP alpha|PITPN

Accession ID NP_006215.1

Blocking Peptide EBP06319

Immunogen Peptide with sequence C-RQKDPVKGMTADD, from the C Terminus of the protein sequence according to NP_006215.1.

Peptide Sequence C-RQKDPVKGMTADD

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat, Dog, Cow

Reactive Species Human

Human Gene ID 5306

Mouse Gene ID 18738

Rat Gene ID 29525

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Limit Antibody detection limit dilution 1:4000.

Western Blot Approx 32kDa band observed in 293 and Human Liver lysates (predicted MW of 31kDa according to NP_006215). Recommended for use at 1-3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

