

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PITPN / PITP ALPHA ANTIBODY

SKU: EB06319



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

phosphatidylinositol transfer protein|phosphotidylinositol transfer protein|PITPNB|PITP alpha|PITPN

Names

Accession

ID

NP_006215.1

Blocking

Peptide

EBP06319

Peptide with sequence C-RQKDPVKGMTADD, from the C Terminus of the protein sequence according to **Immunogen**

NP 006215.1.

Peptide Sequence

C-RQKDPVKGMTADD

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human, Mouse, Rat, Dog, Cow

Reactive

Species

Human

Human

5306

Gene ID Mouse

18738

Gene ID

Rat Gene ID 29525

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade ELISA

Detection

Antibody detection limit dilution 1:4000.

Limit

Western Approx 32kDa band observed in 293 and Human Liver lysates (predicted MW of 31kDa according to

Blot NP_006215). Recommended for use at $1-3\mu g/ml$.

Application

Pep-ELISA, WB

Type





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa