

GOAT ANTI-PJA1 ANTIBODY

SKU: EB05929



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	OTTHUMP00000158770 praja ring finger 1 praja 1 FLJ11830 RNF70 PJA1
Accession ID	NP_660095.1; NP_071763.2
Blocking Peptide	EBP05929
Immunogen	Peptide with sequence QEESKPVWPNT-C, from the N Terminus of the protein sequence according to NP_660095.1; NP_071763.2.
Product Comments	This antibody is expected to recognise both reported isoforms.
Peptide Sequence	QEESKPVWPNT-C
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	64219
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx. 25-30kDa band observed in human brain and human testis lysates at 0.5ug of primary. Please note that the band we observe is different from the predicted size of approx. 74kDa and 52kDa according to AF262024 and NP_071763, respectively. Currently we cannot find an explanation for this in the literature. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any splice variants/modified forms been reported?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)