

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PJA1 ANTIBODY

SKU: EB05929

SPECIFICATIONS

Unit Size 100 µg

Synonym /

Alias OTTHUMP00000158770|praja ring finger 1|praja 1|FL|11830|RNF70|PJA1

Names

ID

Accession

NP_660095.1; NP_071763.2

Blocking

EBP05929 **Peptide**

Peptide with sequence GQESSKPVWPNPT-C, from the N Terminus of the protein sequence according to **Immunogen**

NP 660095.1; NP 071763.2.

Product

This antibody is expected to recognise both reported isoforms. **Comments**

Peptide

GQESSKPVWPNPT-C Sequence

Refrigerated Instructions

Predicted

Human **Species**

Human

64219

Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection

Western

Antibody detection limit dilution 1:32000.

Limit

Blot

Approx. 25-30kDa band observed in human brain and human testis lysates at 0.5ug of primary. Please note that the band we observe is different from the predicted size of approx. 74kDa and 52kDa according to AF262024 and NP_071763, respectively. Currently we cannot find an explanation for this in the literature. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any

splice variants/modified forms been reported?

Application

Pep-ELISA Туре