


GOAT ANTI-PKD1 (AA2281-2292) ANTIBODY

SKU: EB11002



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	polycystin-1 polycystin 1 polycystic kidney disease-associated protein autosomal dominant polycystic kidney disease 1 protein TRPP1 Pc-1 PBP polycystic kidney disease 1 (autosomal dominant) PKD1
Usage Summary	Flow Cytometry: Flow cytometric analysis of HEK293 cells. Recommended concentration: 10ug/ml.
Accession ID	NP_001009944.2; NP_000287.3
Blocking Peptide	EBP11002
Immunogen	Peptide with sequence C-KSYDPNLEDGDQT, from the internal region of the protein sequence according to NP_001009944.2; NP_000287.3.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_001009944.2; NP_000287.3).
Peptide Sequence	C-KSYDPNLEDGDQT

Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig
Reactive Species	Human
Human Gene ID	5310
Mouse Gene ID	18763
Rat Gene ID	24650
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA, FC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

