

GOAT ANTI-PL6 / TMEM115 ANTIBODY

SKU: EB05654



SPECIFICATIONS

Unit Size	100 µg
Synonym /	
Alias	PL6 protein placental protein 6 PL6 transmembrane protein 115 TMEM115
Names	
Accession ID	NP_008955.1
Blocking Peptide	EBP05654
Immunogen	Peptide with sequence C-SSLITFEAAPPTL, from the C Terminus of the protein sequence according to NP_008955.1.
Peptide Sequence	C-SSLITFEAAPPTL
Shipping Instructions	Refrigerated
Predicted Species	Human, Cow
Human Gene ID	11070
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	No signal obtained yet but low background observed in A549, 293 and human kidney lysates at up to 1µg/ml.
Application Type	Pep-ELISA