


GOAT ANTI-PLASMINOGEN ACTIVATOR / PLAU ANTIBODY

SKU: EB06922



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names urokinase-type plasminogen activator|urokinase-type plasminogen activator precursor|urokinase-type plasminogen activator amino-terminal fragment|urokinase plasminogen activator|plasminogen activator, urinary|antagonist of uPA|U-plasminogen activator|u-PA|URK|UPA|ATF|HGNC:9052|plasminogen activator, urokinase|PLAU

Accession ID NP_002649.1; NP_001138503.1

Blocking Peptide EBP06922

Immunogen Peptide with sequence C-HFLPWIRSHTKEEN, from the C Terminus of the protein sequence according to NP_002649.1; NP_001138503.1.

Product Comments This antibody is expected to recognise both reported isoforms.

Peptide Sequence C-HFLPWIRSHTKEEN

Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	5328
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Prostate and Placenta. Recommended concentration: 8µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

