

GOAT ANTI-PLCG2 ANTIBODY

SKU: EB05109



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Phospholipase C, gamma 2 (phosphatidylinositol-specific) phospholipase C, gamma 2 (phosphatidylinositol-specific) PLCG2
Blocking Peptide	EBP05109
Immunogen	Peptide with sequence NHSSCTRRNATRG, from the C Terminus of the protein sequence according to .
Peptide Sequence	NHSSCTRRNATRG
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	5336
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Jurkat and Daudi lysates at upto 2µg/ml.