



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PLCG2 ANTIBODY

SKU: EB05109

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Phospholipase C, gamma 2 (phosphatidylyinositol-specific)|phospholipase C, gamma 2 (phosphatidylinositol-Alias

specific)|PLCG2 Names

Blocking

EBP05109 **Peptide**

Immunogen Peptide with sequence NHSSCTRRNATRG, from the C Terminus of the protein sequence according to .

Peptide

NHSSCTRRNATRG

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Instructions Refrigerated

Predicted

Human **Species**

Human 5336

Gene ID Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

Western

Western blot: No signal obtained yet but low background observed in Jurkat and Daudi lysates at upto 2µg/ml. Blot