Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-PNKP / PNK ANTIBODY**

**SKU:** EB05517

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Homo sapiens polynucleotide kinase 3\1-phosphatase (PNKP)|polynucleotide kinase 3-prime-**Alias** 

phosphatase|PNK|polynucleotide kinase 3\1-phosphatase|PNKP **Names** 

**Accession** NP\_009185.2 ID

**Blocking** EBP05517

**Peptide** 

Peptide with sequence EPRLGRLYCQFSEG, from the C Terminus of the protein sequence according to **Immunogen** 

NP 009185.2.

**Peptide** 

Sequence

**EPRLGRLYCQFSEG** 

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Instructions

Refrigerated

**Predicted** 

**Species** 

Human

Reactive

Human **Species** 

Human

11284

Gene ID **Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:16000.

Limit

Western

Sucessfully used to detect recombinant human PNKP (1/1000 dilution). **Blot** 

Application

**Type** 

Pep-ELISA, WB





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600