


GOAT ANTI-POU4F2 ANTIBODY

SKU: EB10033



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	POU4F2 POU domain, class 4, transcription factor 2 POU domain protein POU domain class 4 transcription factor 2 POU class 4 homeobox 2 Brn3b POU domain transcription factor Brn-3b BRN3B BRN3.2
Usage Summary	Flow Cytometry: Flow cytometric analysis of MCF7 cells. Recommended concentration: 10ug/ml.
Accession ID	NP_004566.2
Blocking Peptide	EBP10033
Immunogen	Peptide with sequence C-EKSHREKLTKE, from the internal region of the protein sequence according to NP_004566.2.
Product Comments	This antibody is not expected to cross-react to the similar POU4F3 (GeneID 5459)

Peptide Sequence	C-EKSHREKLTKE
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	5458
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Application Type	Pep-ELISA, FC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

