



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-POU6F1 ANTIBODY

SKU: EB10040



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias

POU6F1|TCFB1|POU domain, class 6, transcription factor 1|POU class 6 homeobox 1|MPOU|BRN5

Names

ID

Accession

NP_002693.3

Blocking

Peptide

EBP10040

Immunogen

Peptide with sequence C-ETPTVSQLVSKPHTP, from the internal region of the protein sequence according to

NP 002693.3.

Peptide

C-ETPTVSQLVSKPHTP

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method

using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat, Dog, Cow

Reactive

Species

Human, Mouse

5463

Human

Gene ID

Product

Grade

ELISA

Detection

Limit

Antibody detection limit dilution 1:128000.

Western

Blot

Approx 30kDa band observed in Human Brain (Cerebral cortex) and Mouse fetal Brain lysates (calculated MW

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

of 32.6kDa according to NP_002693.3). Recommended concentration: 0.3-1µg/ml.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES





 $\textbf{Email:} \ \underline{\textbf{customerservice@vectorlabs.com}}$

Telephone: (650) 697-3600

250kDa 150kDa		250kDa 150kDa
100kDa		100kDa
75kDa		75kDa
50kDa		50kDa
37kDa		37kDa
25kDa		051-D-
20kDa		25kDa
		20kDa
15kDa		
	1	15kDa