


GOAT ANTI-PPARG ANTIBODY

SKU: EB07241



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / ppar gamma2|peroxisome proliferator-activated receptor gamma 1|peroxisome proliferator activated-receptor

Alias gamma|peroxisome proliferative activated receptor gamma|PPAR

Names gamma|PPARG2|PPARG1|NR1C3|HUMPPARG|gamma|peroxisome proliferative activated receptor|PPARG

Usage Immunofluorescence: Strong expression of the protein seen in A549 and A431 cells.

Summary Recommended concentration: 10µg/ml.

Accession ID NP_005028.5; NP_056953.2; NP_001341597.1; NP_001361194.1; NP_001341599.1; NP_001317544.2; NP_001361195.1

Blocking Peptide EBP07241

Peptide with sequence CDLNCRIHKKSRNK, from the internal region of the protein sequence according to

Immunogen NP_005028.5; NP_056953.2; NP_001341597.1; NP_001361194.1; NP_001341599.1; NP_001317544.2; NP_001361195.1.

Product Comments This antibody is expected to recognise multiple isoforms.

Peptide Sequence CDLNCRIHKKSRNK

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat

Reactive Species Human

Human Gene ID 5468

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:64000.

Application Type Pep-ELISA, IF

GALLERY IMAGES

