


GOAT ANTI-PPIG ANTIBODY

SKU: EB09674



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa


SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PPiG peptidylprolyl isomerase G (cyclophilin G) CARS-Cyp CYP MGC133241 SRCyp Clk-associating RS-cyclophilin peptidyl-prolyl isomerase G (cyclophilin G) peptidylprolyl isomerase G
Accession ID	NP_004783.2
Blocking Peptide	EBP09674
Immunogen	Peptide with sequence C-RDHNSSNSREKKAD, from the internal region of the protein sequence according to NP_004783.2.
Peptide Sequence	C-RDHNSSNSREKKAD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	9360
Mouse Gene ID	228005
Rat Gene ID	83624
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Approx 85kDa band observed in Human Kidney lysates (calculated MW of 88.6kDa according to NP_004783.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa