



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PPIG ANTIBODY

SKU: EB09674



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

PPIG|peptidy|prolyl isomerase G (cyclophilin G)|CARS-Cyp|CYP|MGC133241|SRCyp|Clk-associating RS-**Alias**

cyclophilin|peptidyl-prolyl isomerase G (cyclophilin G)|peptidylprolyl isomerase G Names

Accession

NP 004783.2 ID

Blocking

EBP09674 **Peptide**

Peptide with sequence C-RDHNSSNNSREKKAD, from the internal region of the protein sequence according to **Immunogen**

NP 004783.2.

Peptide C-RDHNSSNNSREKKAD Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog, Cow **Species**

Reactive Human

Species

Human 9360

Gene ID

Mouse 228005 **Gene ID**

Rat Gene ID 83624

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade ELISA

Detection Antibody detection limit dilution 1:1000.

Limit

Western Approx 85kDa band observed in Human Kidney lysates (calculated MW of 88.6kDa according to NP 004783.2).

Blot Recommended concentration: 1-3µg/ml.

Application

Pep-ELISA, WB **Type**

> Goat Anti-PPIG Antibody https://everestbiotech.com/products/goat-anti-ppig-antibody/



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa