


GOAT ANTI-PRDM14 (AA533-544) ANTIBODY

SKU: EB10762



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names PRDM14| PR domain containing 14| MGC59730| PFM11| PR domain zinc finger protein 14| PR domain-containing protein 14| PR-domain zinc finger protein 14

Accession ID NP_078780.1

Blocking Peptide EBP10762

Immunogen Peptide with sequence C-DSHVRRSHKEDD, from the internal region of the protein sequence according to NP_078780.1.

Peptide Sequence C-DSHVRRSHKEDD

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat, Dog, Pig, Cow

Reactive Species Human

Human Gene ID 63978

Mouse Gene ID 383491

Rat Gene ID 100365184

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:2000.

Western Blot Approx 70kDa band observed in lysates of cell line A549 (calculated MW of 64.1kDa according to NP_078780.1). Recommended concentration: 1-3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

