

GOAT ANTI-PRDM5 ANTIBODY

SKU: EB05714



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PFM2 PR domain containing 5 PRDM5
Accession ID	NP_061169.2
Blocking Peptide	EBP05714
Immunogen	Peptide with sequence C-KVHMDNIHGVADS, from the C Terminus of the protein sequence according to NP_061169.2.
Peptide Sequence	C-KVHMDNIHGVADS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Monkey
Human Gene ID	11107
Mouse Gene ID	70779
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Western Blot: No signal obtained yet but low background observed in Human Brain, Human Lungs, and Raji lysates at upto 2µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

