Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PRDM5 ANTIBODY

SKU: EB05714

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias PFM2|PR domain containing 5|PRDM5

Names

Accession NP 061169.2 ID

Blocking

EBP05714 **Peptide**

Peptide with sequence C-KVHMDNIHGVADS, from the C Terminus of the protein sequence according to Immunogen

NP 061169.2.

Peptide C-KVHMDNIHGVADS

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Monkey **Species**

Human

11107 Gene ID

Mouse

70779

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:2000.

Limit

Western Blot

Western Blot: No signal obtained yet but low background observed in Human Brain, Human Lungs, and Raji lysates at upto 2µg/ml. We would appreciate any feedback from people in the field - have any results been reported with

other antibodies/lysates?

Application

Pep-ELISA **Type**





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600