

GOAT ANTI-PRDM8 ANTIBODY

SKU: EB05716



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	PFM5 PR-domain containing protein 8 PR domain containing 8 PRDM8
Accession ID	NP_064611.3; NP_001092873.1
Blocking Peptide	EBP05716
Immunogen	Peptide with sequence EDTGIQRGIWDG-C, from the N Terminus of the protein sequence according to NP_064611.3; NP_001092873.1.
Product Comments	The variants (NP_064611.3; NP_001092873.1) encode the same protein.
Peptide Sequence	EDTGIQRGIWDG-C
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Dog
Human Gene ID	56978
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Preliminary experiments gave no signal but low background in Human Heart and Jurkat lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA