


GOAT ANTI-PRIM1 (RAT) ANTIBODY

SKU: EB10734



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	
Names	DNA primase small subunit, 49kDa MGC109113 DNA primase, p49 subunit Prim1
Accession ID	NP_001008768.1
Blocking Peptide	EBP10734
Immunogen	Peptide with sequence C-RKGALLKKSDLQKD, from the C Terminus of the protein sequence according to NP_001008768.1.
Peptide Sequence	C-RKGALLKKSDLQKD
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Human, Mouse
Mouse Gene ID	19075
Rat Gene ID	246327
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 30kDa band observed in Human and Mouse Liver lysates (calculated MW of 34.6kDa according to NP_001008768.1). Recommended concentration: 0.01-0.03µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

