

GOAT ANTI-PRKAG3 ANTIBODY

SKU: EB08471



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	AMPK gamma-3 chain AMP-activated protein kinase, non-catalytic gamma-3 subunit 5\1-AMP-activated protein kinase, gamma-3 subunit AMPKG3 protein kinase, AMP-activated, gamma 3 non-catalytic subunit PRKAG3
Accession ID	NP_059127.2
Blocking Peptide	EBP08471
Immunogen	Peptide with sequence C-KLGWDDELKPGAQ, from the internal region of the protein sequence according to NP_059127.2.
Peptide Sequence	C-KLGWDDELKPGAQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Cow
Human Gene ID	53632
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)