

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-PRKAG3 ANTIBODY**

**SKU:** EB08471

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

**Alias** 

AMPK gamma-3 chain|AMP-activated protein kinase, non-catalytic gamma-3 subunit|5\1-AMP-activated protein kinase, gamma-3 subunit|AMPKG3|protein kinase, AMP-activated, gamma 3 non-catalytic subunit|PRKAG3

**Names** 

Accession NP 059127.2 ID

**Blocking** 

**Peptide** 

EBP08471

Peptide with sequence C-KLGWDDELRKPGAQ, from the internal region of the protein sequence according to Immunogen

NP 059127.2.

**Peptide** 

Sequence

C-KLGWDDELRKPGAQ

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Instructions

Refrigerated

**Predicted** 

**Species** 

Human, Cow

Human

53632 Gene ID

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

**ELISA** 

**Detection** 

Antibody detection limit dilution 1:8000.

Limit

**Application** Pep-ELISA **Type**