

GOAT ANTI-PRL3 / PTP4A3 ANTIBODY

SKU: EB05302



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	potentially prenylated protein tyrosine phosphatase PRL-R PRL-3 protein tyrosine phosphatase type IVA, member 3 PRL3 PTP4A3
Accession ID	NP_009010.2; NP_116000.1
Blocking Peptide	EBP05302
Immunogen	Peptide with sequence NRPAPVEVSYKHMR-C, from the N Terminus of the protein sequence according to NP_009010.2; NP_116000.1.
Product Comments	This antibody is expected to recognise both reported isoforms (NP_009010.2 and NP_116000.1).
Peptide Sequence	NRPAPVEVSYKHMR-C
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	11156
Mouse Gene ID	19245
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 22kDa band observed in Human Pancreas lysates (calculated MW of 19.5Da according to NP_116000.1). Recommended concentration: 0.03-0.1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES



