

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PRL3 / PTP4A3 ANTIBODY

SKU: EB05302



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

Telephone: (650) 697-3600





SPECIFICATIONS

Unit Size

100 µg

Synonym / Alias

potentially prenylated protein tyrosine phosphatase|PRL-R|PRL-3|protein tyrosine phosphatase type IVA,

member 3|PRL3|PTP4A3 **Names**

Accession

NP_009010.2; NP_116000.1

Blocking

EBP05302 **Peptide**

Immunogen

Peptide with sequence NRPAPVEVSYKHMR-C, from the N Terminus of the protein sequence according to

NP 009010.2; NP 116000.1.

Product Comments

This antibody is expected to recognise both reported isoforms (NP_009010.2 and NP_116000.1).

Peptide

NRPAPVEVSYKHMR-C

Sequence

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat

Reactive

Species

Human

Human

11156 **Gene ID**

Mouse

19245 **Gene ID**

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Antibody detection limit dilution 1:64000.

Limit

Approx 22kDa band observed in Human Pancreas lysates (calculated MW of 19.5Da according to Western

Blot NP_116000.1). Recommended concentration: 0.03-0.1µg/ml.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa