

GOAT ANTI-PRODH (AA170-184) ANTIBODY

SKU: EB09081



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	tumor protein p53 inducible protein 6 proline oxidase, mitochondrial proline oxidase 2 proline dehydrogenase 1 proline dehydrogenase (proline oxidase) p53 induced protein TP53I6 SCZD4 PRODH2 PRODH1 POX PIG6 MGC148079 MGC148078 HSPOX2 FLJ33744 proline dehydrogenase (oxidase) 1 PRODH
Accession ID	NP_057419.3
Blocking Peptide	EBP09081
Immunogen	Peptide with sequence C-DGSGTNKRDKQYQAH, from the internal region of the protein sequence according to NP_057419.3.
Peptide Sequence	C-DGSGTNKRDKQYQAH
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	5625
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Preliminary experiments in Human Brain (Cerebellum, Frontal Cortex and Substantia Nigra) lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA