

GOAT ANTI-PRPF6 ANTIBODY

SKU: EB09694



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PRPF6 PRP6 pre-mRNA processing factor 6 homolog (S. cerevisiae) ANT-1 C20orf14 TOM U5-102K hPrp6 OTTHUMP00000031625 PRP6 pre-mRNA processing factor 6 homolog U5 snRNP-associated 102 kDa protein androgen receptor N-terminal domain transactivating protein-1 p102 U5 small nuclear ribonucleoprotein particle-binding protein putative mitochondrial outer membrane protein import receptor
Accession ID	NP_036601.2
Blocking Peptide	EBP09694
Immunogen	Peptide with sequence C-TGENHTSVDPRQTQ, from the Internal region of the protein sequence according to NP_036601.2.
Peptide Sequence	C-TGENHTSVDPRQTQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	24148
Mouse Gene ID	68879
Rat Gene ID	366276
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 110kDa band observed in nuclear lysates of cell line Jurkat (calculated MW of 107kDa according to NP_036601.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

