

GOAT ANTI-PSD / TYL ANTIBODY

SKU: EB06172



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	PSD1 KIAA2011 RP11-18I14.2 pleckstrin and Sec7 domain containing TYL pleckstrin and Sec7 domain protein PSD
Accession ID	NP_002770.1
Blocking Peptide	EBP06172
Immunogen	Peptide with sequence PLKSPVPFLPGTS-C, from the N Terminus of the protein sequence according to NP_002770.1.
Peptide Sequence	PLKSPVPFLPGTS-C
Shipping Instructions	Refrigerated
Predicted Species	Human, Rat, Dog
Human Gene ID	5662
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Preliminary experiments gave no signal but low background in Human Brain extracts at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA