

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-PSPHL / CO9 ANTIBODY**

**SKU:** EB06114

×

## **SPECIFICATIONS**

Unit Size 100 μg

Synonym /

Alias L-3-phosphoserine-phosphatase homolog|CO9|phosphoserine phosphatase-like|PSPHL

Names

Accession

NP\_003823.1 removed

Blocking Peptide

EBP06114

Peptide with sequence HWIWRKCDQATSQG, from the C Terminus of the protein sequence according to

NP 003823.1 removed.

Peptide

Sequence HWIWRKCDQATSQG

Shipping

Instructions Refrigerated

Predicted

Human

Species

Human 8781

Ochic ib

Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:16000.

Limit

Lilling

Western Blot Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the predicted size of approx. 7.80 kDa according to NP\_003823. Therefore we cannot recommend an optimal

concentration and the product is investigative grade.

**Application** 

Type Pep-ELISA