

GOAT ANTI-PSPHL / CO9 ANTIBODY

SKU: EB06114



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	L-3-phosphoserine-phosphatase homolog CO9 phosphoserine phosphatase-like PSPHL
Accession ID	NP_003823.1 removed
Blocking Peptide	EBP06114
Immunogen	Peptide with sequence HWIWRKCDQATSQG, from the C Terminus of the protein sequence according to NP_003823.1 removed.
Peptide Sequence	HWIWRKCDQATSQG
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	8781
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the predicted size of approx. 7.80kDa according to NP_003823. Therefore we cannot recommend an optimal concentration and the product is investigative grade.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

