




Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PTCD2 ANTIBODY

SKU: EB11896







250kDa 150kDa

100kDa

75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Aliquot and store at -20°C. Minimize freezing and thawing.

Accession

NP 079030.3 ID

Blocking

Peptide

EBP11896

Peptide with sequence C-TLHITGQVTTDSLD, from the internal region (near C terminus) of the protein sequence **Immunogen**

according to NP_079030.3.

Peptide

Sequence

C-TLHITGQVTTDSLD

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Human

79810

Predicted

Species

Reactive

Human **Species**

Human

Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Antibody detection limit dilution 1:128000.

Limit

Western

Approx 45kDa band observed in lysates of cell line Daudi (calculated MW of 44.0kDa according to

NP_079030.3). Recommended concentration: 1-3µg/ml. Blot

Application

Type

Pep-ELISA, WB

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa