

GOAT ANTI-PTCHD3 ANTIBODY

SKU: EB10593



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PTCHD3 patched domain containing 3 FLJ44037 MGC129888 PTR patched domain-containing protein 3 patched-related protein
Accession ID	NP_001030014.2
Blocking Peptide	EBP10593
Immunogen	Peptide with sequence C-KKVDYWDKDVRQ, from the internal region of the protein sequence according to NP_001030014.2.
Peptide Sequence	C-KKVDYWDKDVRQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	374308
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 90kDa band observed in Human Testis lysates (calculated MW of 86.9kDa according to NP_001030014.2). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

