

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PTF1A (ZEBRAFISH) ANTIBODY

SKU: EB07973

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias ptf1a|pancreas specific transcription factor, 1a|MGC112216|PTF1|si:zc142h2.2|zgc:112216

Names

Accession NP 997524.1 ID

Blocking

EBP07973 **Peptide**

Peptide with sequence C-EDPRKLHLKSSINN, from the internal region of the protein sequence according to Immunogen

NP 997524.1.

Peptide

C-EDPRKLHLKSSINN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Zebrafish

antibodies/lysates?

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection

Antibody detection limit dilution 1:16000.

Limit

Western

Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We

Blot would appreciate any feedback from people in the field - have any results been reported with other

Application

Pep-ELISA

Type