


GOAT ANTI-PTIP (N TERMINUS) ANTIBODY

SKU: EB10703



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	protein encode PAX transcription activation domain interacting protein 1 like PAX interacting protein 1 TNRC2 PTIP PAXIP1L PACIP1 FLJ41049 CAGF29 CAGF28 PAX interacting (with transcription-activation domain) protein 1 PAXIP1
Names	
Accession ID	NP_031375.3
Blocking Peptide	EBP10703
Immunogen	Peptide with sequence SDQAPKVPEEMFRE-C, from the N Terminus of the protein sequence according to NP_031375.3.
Peptide Sequence	SDQAPKVPEEMFRE-C
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	22976
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 95kDa band observed in lysates of cell line HeLa (calculated MW of 121kDa according to NP_031375.3). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

