

GOAT ANTI-PTPH1 / PTPN3 ANTIBODY

SKU: EB06241



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	cytoskeletal-associated protein tyrosine phosphatase protein tyrosin phosphatase H1 protein tyrosine phosphatase, non-receptor type 3 PTPH1 PTPN3
Accession ID	NP_002820
Immunogen	Peptide with sequence C-RVYEEGLVQMLDPS, from the C Terminus of the protein sequence according to NP_002820.
Peptide Sequence	C-RVYEEGLVQMLDPS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	5774
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Western Blot: Preliminary experiments gave no signal but low background in Hela and A431 lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA