

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PTPH1 / PTPN3 ANTIBODY

SKU: EB06241

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

cytoskeletal-associated protein tyrosine phosphatase|protein tyrosin phosphatase H1|protein tyrosine **Alias**

phosphatase, non-receptor type 3|PTPH1|PTPN3 **Names**

Accession

NP 002820

Peptide with sequence C-RVYEEGLVQMLDPS, from the C Terminus of the protein sequence according to Immunogen NP_002820.

Peptide C-RVYEEGLVQMLDPS Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human 5774

Gene ID

Product Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/aspiring_medium.png$

ELISA

Detection Antibody detection limit dilution 1:8000.

Limit

Western Blot: Preliminary experiments gave no signal but low background in Hela and A431 lysates at up to Western 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other **Blot**

antibodies/lysates?

Application

Pep-ELISA **Type**