

GOAT ANTI-PTPRH ANTIBODY

SKU: EB08377



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PTPRH protein tyrosine phosphatase, receptor type, H FLJ39938, MGC133058, MGC133059, SAP-1
Accession ID	NP_002833.2
Blocking Peptide	EBP08377
Immunogen	Peptide with sequence DHVRKNERDSNC, from the interal region of the protein sequence according to NP_002833.2.
Peptide Sequence	DHVRKNERDSNC
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	5794
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Preliminary experiments in Human Brain (Cerebellum, Cerebral Cortex, Frontal Cortex and Hippocampus) lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)