



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-PTPRH ANTIBODY**

**SKU:** EB08377

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias PTPRH|protein tyrosine phosphatase, receptor type, H|FL|39938, MGC133058, MGC133059, SAP-1

**Names** 

Accession NP 002833.2

**Blocking** 

EBP08377 Peptide

Peptide with sequence DHVRKNERDSNC, from the interal region of the protein sequence according to Immunogen

NP 002833.2.

**Peptide** 

**DHVRKNERDSNC** Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human **Species** 

Human

5794 Gene ID

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:32000.

Limit

Preliminary experiments in Human Brain (Cerebellum, Cerebral Cortex, Frontal Cortex and Hippocampus) lysates Western gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any

Blot feedback from people in the field - have any results been reported with other antibodies/lysates?

Application

Pep-ELISA **Type**