


GOAT ANTI-RAB11FIP1 / RCP ANTIBODY

SKU: EB06897



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	RAB11FIP1 RAB11 family interacting protein 1 (class I) HGNC:30265 FLJ22524 FLJ22622 MGC78448 NOEL1A RCP rab11-FIP1 Rab coupling protein Rab effector protein Rab-interacting recycling protein Rab11-family interacting protein 1D no Rab11-binding protein 1A DKFZp686E2214 RAB11 coupling protein RAB11 family interacting protein 1
Accession ID	NP_079427.3; NP_001002814.1
Blocking Peptide	EBP06897
Immunogen	Peptide with sequence C-EATKEAKESKKPE, from the internal region of the protein sequence according to NP_079427.3; NP_001002814.1.
Product Comments	This antibody is expected to recognise isoform 1 (NP_079427.3) and isoform 3 (NP_001002814.1).
Peptide Sequence	C-EATKEAKESKKPE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	80223
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 80kDa band observed in lysates of cell line HeLa (calculated MW of 71.0kDa according to NP_079427.3). The observed molecular weight of 80kDa corresponds to earlier findings in literature with different antibodies (Lindsay et al, J Biol Chem. 2002 Apr 5;277(14):12190-9. PMID 11786538.). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

