

GOAT ANTI-RAB2 / RAB2A ANTIBODY

SKU: EB05030

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	
Names	GTP binding protein RAB2A RAB2, member RAS oncogene family RAB2A, member RAS oncogene family RAB2
Accession ID	NP_002856.1
Blocking Peptide	EBP05030
Immunogen	Peptide with sequence HAGNQGGQQAGGGCC, from the C Terminus of the protein sequence according to NP_002856.1.
Peptide Sequence	HAGNQGGQQAGGGCC
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Pig, Cow
Reactive Species	Human
Human Gene ID	5862
Mouse Gene ID	59021
Rat Gene ID	65158
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Approx 23kDa and 18kDa bands seen in HeLa cell lysates. Recommended for use at 1-3µg/ml.
Application Type	WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

