

GOAT ANTI-RAB25 ANTIBODY

SKU: EB05498



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	RAB11C CATX-8 RAB25, member RAS oncogene family RAB25
Names	
Accession ID	NP_065120.1
Blocking Peptide	EBP05498
Immunogen	Peptide with sequence QISDCSLFTAPSGS, from the C Terminus of the protein sequence according to NP_065120.1.
Peptide Sequence	QISDCSLFTAPSGS
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	57111
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western blot: Approx 30kDa band seen in Human Kidney lysate (predicted MW of 28kDa according to NP_065120). Recommended for use at 4-6µg/ml.